

New Clinical Evidence for Endoscope processing is in. Implement ANSI/AAMI ST91 Today!

university.STERIS.com

Disclosures

- **Successful completion:** Participants must complete the entire program and submit required documentation. No partial credit will be given.
- **Conflict of interest:** Employee of STERIS.
- **Commercial company support:** Fees are underwritten by education funding provided by STERIS.
- **Non-commercial company support:** None.
- **Alternative/Complementary therapy:** None.

Continuing Education

- STERIS Corporation is an approved provider of continuing education by California Board of Nursing (provider # CEP 11681) and the Board of Ambulatory Surgery Certification (provider # 1417).
- This program is approved for:
 - 1 hour(s) of GI specific content credit by American Board of Certification for Gastroenterology Nurses (ABCGN)
 - 1 AEU(s) & 1 IPCH(s) by BASC
 - 1 contact hour(s) of continuing education credit
 - HSPA • CBSPD
 - CBRN

Continuing Education

The Leader in Perioperative Certification

Through a partnership with CCI®, it also meets CNOR® and CSSM® recertification requirements for perioperative nurses.

Learning Objectives

- Identify key clinical evidence driving the changes in ANSI/AAMI ST91
- Implement ANSI/AAMI ST91 guidelines to develop a comprehensive department for GI, bronchoscopy and other flexible endoscope processing
- Prioritize critical processes and applications for patient safety

FDA Post Market Surveillance

- 4.1- 6.6% contamination rate of duodenoscopes with “high concern microorganisms” (FDA 522 Post Market Surveillance Studies August 2019)
- Gram-negative rods
- Staphylococcus aureus
- Beta-hemolytic Streptococcus species
- Enterococcus species
- Yeasts

Processing Area Design

"Until such time (major renovations or new construction) that two separate processing rooms can be provided, strict unidirectional processing procedures should be in place to reduce risks of cross-contamination."

ANSI/AAMI ST91 Section 4.2.2 pg. 13

"Dirty" Side Design Considerations

Room

- Traffic Control / Doors
- Space
- ANSI/ASHRAE/ASHE 170
- Negative Room Pressure
- Pass thru windows
- Cleanable

Work Area

- Utility water
- Instrument air
- Lighting
- Inspection/Leak test areas
- Hand Washing
- Emergency Eye Wash / Showers

"Clean" Side Design Considerations


Room

- Traffic Control / Doors
- Space
- ANSI/ASHRAE/ASHE 170
- Positive Room Pressure
- Pass thru windows
- Cleanable

Work Area

- Dedicated rinse sink
- Instrument air/HEPA filtered
- Lighting
- Inspection
- Utilities
- Drying
- Hand Washing/ Emergency Eye Wash / Showers


Personnel




The Exposure Control Plan

Standard Precautions	Hand Hygiene	Attire	Personal Protective Equipment (PPE)
• Centers for Disease Control and Prevention (CDC) • OSHA Blood-borne pathogen regulation	• Short fingernails • No artificial nails • CDC and WHO guidelines	• Uniforms provided by the facility • Head and facial hair covered • Non-skid shoes • No jewelry or wristwatches	• Attire • General-purpose utility or chemical resistant gloves • Liquid-resistant coverings, shoe covers and face mask • Eye protection


Point Of Use Treatment



Point Of Use Treatment




Point Of Use Treatment



Transport

- Endoscope and its components
- Closed container or cart
- Labeled Biohazard
- Kept Moist



Handoff Communication



Patient Identifier



Date of Procedure



Completion Time
of Point-of-Use

How Much Time Is Allowed?



"If no timeframe is given, manual cleaning should be initiated within **one hour** or as determined by the facility based on their own documented risk assessment."


ANSI/AAMI ST91, 7.5.1 pg. 36

Delayed Processing


- Cleaning starts with submersion in cleaning solution
- Document time period
- Follow Manufacturer IFU for delayed processing
- Documented procedures

Leak Test

- Detects damage that could expose internal components



Sinks



Manual Cleaning – Prepare

- Fresh compatible cleaning solution
- Video cap attached
- Fresh non-linting cloth or appropriate sponge
- Brushes as defined by endoscope and accessories IFU
- Compatible suction and flushing mechanisms

Manual Cleaning – Clean

- Beneath the cleaning solution surface
- Soft lint-free cloth or appropriate sponge
- Use appropriate brush for ports, openings and elevator housing
- Brush accessible channels
- Aspirate and flush channels per IFU
- Soak as directed



Elevator Guide



Manual Cleaning - Rinse

- Volume of water
- Number of rinses
- Water pressure
- Utility Water per ANSI/AAMI ST108

Common Water Conditions	Standard Specification
Hardness	< 150 Mg/L
Conductivity (mg/L = ppm)	< 500 μ S/cm
pH	6.5-9.5
Total Alkalinity	<400 mg CaCO ₃ /L
Bacteria	<500 CFU/ml
Endotoxin	n/a


*Refer to ST108 for the complete specification

Drying After Cleaning

- Residual could dilute in HLD or LCS process
- Residual could interfere with terminal sterilization
- Designated drying area
- Clean non-linting cloth
- Forced instrument or HEPA filtered air

Lighted Magnification Inspection

- Lighted magnification 5x-10x
- Inspect for:
 - Cleanliness
 - Integrity
 - Moisture
 - Damage



Borescope Inspection

- Inspect accessible channels and openings
- Inspect for:
 - Residual Soils
 - Surface Damage
 - Residual Moisture
 - Manufacturer's inspection criteria
- Competency Program




Cleaning Verification Testing

- Residual Soils - Organic or inorganic
- Not visible
- After cleaning prior to disinfection/sterilization
- Every high-risk endoscope
- Periodically for low-risk endoscopes



Cleaning Verification Program


1. Create standard work, policies & procedures
2. Identify soil markers, benchmarks, and endoscopes
3. Training and competency
4. Perform user verification
5. Implement and report on results



High-Risk Endoscope

"Endoscopes that have been associated with infectious outbreaks including those that are difficult to process and increase the risk of incomplete clearance of contaminating infectious organisms, including bronchoscopes, cystoscopes, duodenoscopes, endobronchial ultrasound endoscopes, linear ultrasound endoscopes, ureteroscopes, and others as determined by the facility."

ANSI/AAMI ST91 3.31 pg. 6



Choosing The Right Biocidal Process

Critical Device	Semi-Critical Device
- Introduced into the blood - Contact sterile tissue or body space - Sterilization	- Contact intact mucous membranes or non-intact skin - Disinfection, preferably sterilization

*"It is advised that flexible and semi-rigid endoscopes to be used in semi-critical applications be sterilized prior to use."
ANSI/AAMI ST91 pg. xiii*

Semi-Critical Endoscope Sterilization

- Based on Risk Assessment
 - High-risk endoscope
 - Patient population
 - Infection rates
 - Outbreak organisms

Manual High-Level Disinfection / Liquid Chemical Sterilization

```

graph TD
    A[Don PPE] --> B[Prepare HLD/LCS Solution]
    B --> C[Submerge and flush lumens]
    C --> D[Remove and purge]
    D --> E[Utility/Critical Water Rinse]
          
```

"Manual processing is not recommended but if necessary due to facility or resource limitations"
ST91 section 8.2.4.1 pg.44

Automated High-Level Disinfection And Liquid Chemical Sterilization

```

graph TD
    A[Don PPE] --> B[Load and attach connectors]
    B --> C[Run Cycle]
    C --> D[Purging]
          
```

"If an endoscope is left in the basin of an AER or LCS for an extended period of time (i.e., more than one hour), the endoscopes should be put back through the HLD/LCS process before it is stored or used"
ST91 section 8.2.3.2 pg.43

Important Considerations For Manual OR Automated HLD And LCS



- Follow IFU
 - Safety
 - Preparation, contact time and rinsing
 - Disposal



- Single vs. Multi Use
 - Maximum reuse
 - No topping-off
 - Test strips for concentration

Important Considerations For Manual OR Automated HLD And LCS



- Chemistry Contact
 - Remove trapped air
 - Appropriate connectors
 - Validated for processor



- Rinse Water Quality
 - Correct number of rinses
 - Manual: Critical (ST108) or Sterile Water
 - Processor: Follow IFU
 - Fresh water for every rinse (by immersion)

When To Dry HLD/LCS Endoscopes

- Prevents microbial growth and biofilm formation during storage
- Semi-critical endoscopes
 - Prior to patient use
- Critical Endoscopes Processed in LCS
 - Use immediately with no drying or dry prior to storage
 - If HLD AER Dry before storage and use

Important Drying Considerations

- AER purge but do not dry
- Use of 70 – 90% alcohol

If the use of alcohol is a consideration, a multidisciplinary team that includes infection preventionists, endoscopy and perioperative RNs, endoscopy processing personnel, endoscopists, and other involved personnel should conduct a risk assessment to determine whether endoscope lumens should be flushed with 70 % to 90 % ethyl or isopropyl alcohol.

ANSI/AAMI ST91 section 8.2.5.1 pg. 46

Endoscope Drying Channels After Cleaning



Minimum

- 10 Minutes
- Drying verification – Annex K

Pressure-Regulated

- Forced Instrument or HEPA Filtered Air
- Regulated to protect endoscope

Dedicated / Sectioned Drying Cabinet

- HEPA-filtered air

Storage

- Closed Cabinet
- HEPA Filtered Air
- Directed HEPA filter air (preferred)
- Storage Time / "hang time"



Terminal Sterilization

Prepare endoscope

Package in sterile barrier

Run Cycle

Release the load



Terminal Sterilization Options

Ethylene Oxide Sterilization

- Validated
- Venting cap and/or removal of cleaning cap
- ANSI/AAMI ST41

Hydrogen Peroxide and Hydrogen Peroxide - Ozone

- Validated
- Venting cap and/or removal of cleaning cap
- ANSI/AAMI ST58



Don't Forget The Accessories!



Quality Control

- Policies and Procedures
- Endoscope and accessories identification and traceability
- Process monitoring and testing
- Preventative maintenance and repairs
- Documentation
- Quality process improvement

Prioritize And Plan Implementation



1. Obtain ANSI/AAMI ST91
2. Audit your process
3. Perform a risk assessment
4. Prioritize - Pareto analysis
5. Make a plan

Action Items

- Obtain copies of the applicable standards and guidelines for flexible endoscope processing for your department
- Create an assessment tool based on ANSI/AAMI ST91
- Audit, prioritize and plan your compliance path to ANSI/AAMI ST91

References

- Association for the Advancement of Medical Instrumentation. (2018). AAMI/ANSI ST58:2013/(R)2018 Chemical sterilization and high-level disinfection in health care facilities. Arlington, VA: AAMI.
- Association for the Advancement of Medical Instrumentation. (2018). AAMI/ANSI ST41:2008/(R)2018 Ethylene oxide sterilization in healthcare facilities: Safety and effectiveness. Arlington, VA: AAMI.
- Association for the Advancement of Medical Instrumentation. (2021). AAMI/ANSI ST91:2021 Flexible and semi-rigid endoscope processing in health care facilities. Arlington, VA: AAMI.
- Association for the Advancement of Medical Instrumentation (2023) AAMI ST108:2023 Water for the processing of medical devices. Arlington, VA: AAMI

References

- Center for Disease Control and Prevention. (2015). *CDC/FDA Health Update about the Immediate Need for Healthcare Facilities to Review Procedures for Cleaning, Disinfecting, and Sterilizing Reusable Medical Devices*. Retrieved from <https://stacks.cdc.gov/view/cdc/34620>
- Kovaleva, J., Peters, F. T. M., van der mai, H.C., and Degener, J.E. (2013). *Transmission of infection by flexible gastrointestinal endoscopy and bronchoscopy*. Clin Micro Rev 26: 231-254. DOI: 10.1128/CMR.00085-12
- Occupational Safety and Health Administration. *1910.1030 Bloodborne Pathogens*. Retrieved from <https://www.osha.gov/laws-regs/regulations/standardnumber/1910/1910.1030>

References

- Ofstead C.L., Heymann OL, Quick MR, Eiland JE, Wetzler HP. Residual moisture and waterborne pathogens inside flexible endoscopes: evidence from a multisite study of endoscope drying effectiveness. Am J Infect Control. 2018;46:689-696.
- Ofstead, C. L., Dirlam Langlay, A. M., Mueller, N. J., Tosh, P. K., & Wetzler, H. P. (2013). Re-evaluating endoscopy-associated infection risk estimates and their implications. American Journal of Infection Control, 41(8), 734-736.
- Wang P, Xu T, Ngamruengphong S, Makary MA, Kallou A, Hutfless S. Rates of infection after colonoscopy and esophagogastroduodenoscopy in ambulatory surgery centres in the USA. Gut. 2018 Sep;67(9):1626-1636.
- Society of Gastroenterology Nurses and Associates, I. (2016). Standards of infection prevention in reprocessing flexible gastrointestinal endoscopes. Chicago, IL: SGNA.



Questions?



Thank you!